

# Regenix™ Liver

## PROTOCOL

**Regenix™ Liver** is composed of various basement membrane proteins separated from liver tissues. Regenix™ Liver can be utilized for two-dimensional (2D) and three-dimensional (3D) culture of hepatic cells. In particular, Regenix™ Liver can provide an optimized environment for adult stem cell (AdSC)-derived and pluripotent stem cell (PSC)-derived liver organoids.

### PROCEDURE

#### 3D culture of liver organoid using Regenix™ Liver

- 01

Thaw Regenix™ Liver for at least 4 hours by submerging the vial in an ice bucket and storing it in a 4°C refrigerator before use. Avoid multiple freeze/thaw cycles.
- 02

Mix Regenix™ Liver by slowly pipetting; Be careful not to create air bubbles during this process.

**Note** Regenix™ Liver may form a gel at temperature above 10°C. The temperature must be lowered to between 4°C and 8°C throughout all handling processes to ensure depolymerization.
- 03

Add Regenix™ Liver to the cell pellet and resuspend evenly by slow pipetting.

**Note** It is recommended to remove as much of the supernatant as possible before adding the Regenix™ Liver.
- 04

Dispense 30 µL of the mixture to each well of a 48-well plate, and then incubate at 37 °C for 40 mins.
- 05

Add the appropriate volume of medium very slowly.

**Note** If you need to add 300 µL per well medium to each well, add the medium slowly and carefully over 15 seconds.

**Note** Culture of liver organoids with Regenix™ Liver requires the addition of 10µM Y-27632 in the first 1-2 days.

# Regenix™ Liver

## PROTOCOL

### PROCEDURE

#### Passage of liver organoids in Regenix™ Liver

- 01** • Prepare 2 mg/mL of collagenase IV (600 – 800 U/mL) in basal medium.  
**Note** Different types of collagenase also work, but need to be optimized to the proper concentration.
- 02** • Gently touch the side of Regenix™ Liver droplet with a 1000 µL pipette tip to detach it from the bottom of the well plate.
- 03** • Cut the tip off a 1000 µL pipette tip with sterile scissors to obtain an opening of 2.5-3 mm in diameter, and use it to transfer each Regenix™ Liver encapsulating organoids to a 15 mL conical tube.  
**Note** It is recommended to use a 15 mL conical tube to avoid cell pellet sticking to the microtube wall.
- 04** • Gently aspirate the supernatant and add enough collagenase IV solution to fully submerge Regenix™ Liver droplets. (e.g. Use 1 mL collagenase IV solution per 6-8 Regenix™ Liver droplets.
- 05** • Incubate the 15 mL conical tube containing Regenix™ Liver upright in a 37 °C incubator for 1 hour.  
**Note** Incubation times longer than 1 hour can damage the organoids.
- 06** • After 1 hour, a thin layer of Regenix™ Liver above the organoid pellet can be seen. Carefully aspirate the layer of Regenix™ Liver and wash the organoids twice with basal medium.
- 07** • Re-encapsulate the organoids in Regenix™ Liver and cultivate them in the same way as before.

### CAUTION

\* Organoids cultured in Regenix™ Liver grow slower than in Matrigel, but there is no observable difference in final organoid size. It is recommended to set the analysis/ harvest timepoint 2-3 days later than when using Matrigel. (e.g. Organoids grown in Regenix™ Liver for 7 days have the same observable size as organoids grown in Matrigel for 5 days )