

# Regenix™ Xenograft

## PROTOCOL

**Regenix™ Xenograft** is designed for patient-derived (PDX) and cell line-derived (CDX) xenograft applications. Composed of various basement membrane proteins from normal porcine tissues and featuring atelocollagen for biological safety, it offers easy sample preparation and injection. Regenix™ Xenograft ensures high cell retention, promoting rapid tumor growth with a smooth, round shape, making it a reliable tool for consistent, high-quality research results.

### PROCEDURE

#### CDX model generation using Regenix™ Xenograft

- 01** Thaw Regenix™ Xenograft for at least 4 hours by submerging the vial in an ice bucket and storing it in a 4°C refrigerator before use. Avoid multiple freeze/thaw cycles.
- 02** Cut the tip off a 200 µL pipette tip with sterile scissors to obtain an opening diameter of 1.5-2 mm, and mix Regenix™ Xenograft by slowly pipetting; Be careful not to create air bubbles during this process.
  - Note** Regenix™ Xenograft may have high viscosity, so there may be some difficulty in pipetting. If a lot of bubbles are generated after pipetting, centrifuge before use.
  - Note** Regenix™ Xenograft may form a gel at temperature above 10°C. The temperature must be lowered to between 4°C and 8°C throughout all handling processes to ensure depolymerization.
- 03** Add Regenix™ Xenograft to the cell pellet and resuspend evenly by slow pipetting.
  - Note** Keep the mixture cold to prevent the Regenix™ Xenograft from solidifying before injection.
- 04** Using a chilled syringe with a 25-29 gauge needle, draw up the 100 µL cell-Regenix™ Xenograft mixture.
- 05** Inject the cell-Regenix™ Xenograft mixture subcutaneously into the flank of the anesthetized mouse. Slowly administer the entire volume of the mixture.
  - Note** Insert the needle at a shallow angle to ensure the injection remains in the subcutaneous space and does not penetrate the muscle.