

# Regenix Pancreas protocol



Regenix Pancreas is composed of various basement membrane proteins separated from the pancreas tissues. Regenix Pancreas can be utilized for two-dimensional (2D) and three-dimensional (3D) culture of pancreatic cells. In particular, Regenix Pancreas can provide an optimized environment for adult stem cells (AdSCs)-derived pancreatic organoids.

## Storage Instructions

- ⊙ Avoid storing Regenix Pancreas on freezer doors or in frequently opened freezers.
- ⊙ After the initial thaw, aliquot Regenix Pancreas into freezer-compatible tubes and store at -80°C. Minimize repeated freezing and thawing to maintain product quality.
- ⊙ Long-term storage after thawing is not recommended for optimal product integrity.
- ⊙ Frozen Regenix Pancreas is stable for up to 2 years from the date of manufacture.

## Thawing Instructions

- ⊙ Regenix Pancreas begins to gel at temperatures above 10°C.
- ⊙ Thaw for at least 4 hours at 2°C to 8°C, ensuring the vial is fully surrounded by ice.
- ⊙ During thawing, keep the ice bucket covered and place it in a cold room or at the back of a refrigerator for consistent temperature control.

## Instructions for 3D Culture of Pancreatic Organoids

### ⊙ Preparation of Regenix Pancreas

Thaw Regenix Pancreas and gently mix by slow pipetting. Due to its high viscosity, use a 200 µL pipette tip cut to a 1.5–2 mm opening to reduce bubble formation. If bubbles occur, centrifuge before use. Keep Regenix Pancreas at 4–8°C during handling to prevent gelation above 10°C.

### ⊙ Organoid Resuspension

Before adding Regenix Pancreas, carefully remove as much supernatant as possible from the prepared organoid pellet. Then, add Regenix Pancreas and gently mix by slow pipetting to ensure uniform resuspension. Regenix Pancreas is provided as a ready-to-use pre-gel solution. Dilution is not recommended, as it may prevent proper hydrogel formation.

### ⊙ Gelation

Dispense 30 µL of the mixture into each well of a 48-well plate and incubate at 37°C for 40 minutes to allow gel formation.

## ☉ Medium Addition

Carefully add the appropriate volume of medium. If adding 300  $\mu$ L per well, dispense the medium slowly over 15 seconds to prevent disruption of the gel. The culture of organoids with Regenix Pancreas may require the addition of 10  $\mu$ M Y-27632 during the first 1–2 days.

## Subculture Instructions for Pancreatic Organoids

### ☉ Preparation of Collagenase Solution

Prepare a 2 mg/mL solution of collagenase IV (600–800 U/mL) in basal medium. Other types of collagenase can be used, but their concentrations may need optimization.

### ☉ Detachment of Regenix Pancreas Droplets

Gently touch the side of the Regenix Pancreas droplet with a 1000  $\mu$ L pipette tip to detach it from the bottom of the well plate.

### ☉ Transfer of Regenix Pancreas-Encapsulated Organoids

Cut off the tip of a 1000  $\mu$ L pipette tip with sterile scissors to create a 2.5–3 mm opening. Use this modified tip to transfer the Regenix Pancreas-encapsulated organoids into a 15 mL conical tube. Using a 15 mL conical tube is recommended to prevent the organoid pellet from sticking to the microtube walls.

### ☉ Collagenase Treatment

Carefully remove the supernatant and add enough collagenase IV solution to fully submerge the Regenix Pancreas droplets (e.g., use 1 mL of collagenase IV solution per 6–8 Regenix Pancreas droplets).

Incubate the 15 mL conical tube upright at 37°C for 1 hour. Do not exceed 1 hour, as prolonged incubation may damage the organoids.

### ☉ Removal of Regenix Pancreas and Washing

After 1 hour, a thin layer of Regenix Pancreas will remain above the organoid pellet. Carefully aspirate this layer and wash the organoids twice with basal medium.

### ☉ Re-Encapsulation

Re-encapsulate the organoids in Regenix and continue culturing under the same conditions as before.